03

Serie SK07

AUTOCLAVE (16 lt - 21 lt)

AUTOCLAVE (16 lt - 21 lt)

MANUALE D USO

INSTRUCTION MANUAL

MODE D EMPLOI

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MANUAL DE USO

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CONTENTS

Safety standards	page 34
Characteristics	page 36
Description of parts*	page 38
Installation.....	page 40
Instructions for use.....	page 42
Maintenance.....	page 54
Guide to problems.....	page 60
Technical specifications	page 62

* The images are reported on the last pages

SYMBOLS

DANGER

The paragraphs marked with this symbol contain instructions which must be carefully carried out to avoid damaging the device, harming the operator and possibly the patient.

WARNINGS

These instructions warn you to pay utmost attention to avoid situations that could damage the device.

PROHIBITION

This icon highlights what should be avoided to prevent damaging the device.

SUGGESTIONS

This icon provides information that allows you to use the device more efficiently.

Disposal container. Dispose the device observing the standards about the differentiated collection of electrical material.

WARNING

(yellow background)

Comply with the instructions outlined in this manual.

WARNING

Symbol of risk or danger. Refer to the manual.

WARNING

Hot parts.

PACKAGING SYMBOLS

HIGH

FRAGILE

DO NOT WET


DO NOT ROLL

DO NOT USE HOOKS

MAX SUPPORTED WEIGHT

**STORAGE TEMPERATURE
TRANSPORT CONDITION**

RH
10% ÷ 90%

RELATIVE HUMIDITY

P
500 hPa ÷ 1060 hPa

ATMOSPHERIC PRESSURE

RECYCLABLE CARDBOARD

SAFETY STANDARDS

Follow all the instructions outlined in this manual.

- The autoclave must be used by staff specifically trained in the sterilisation process, in medical studios, chiropodist studios and beauty salons under the responsibility of the Director of Health.

- Device intended to sterilise, with steam at a temperature of 121°C or 134°C, tools, rotating instruments, re-usable medical devices and porous materials (e.g. fabric).

- Only insert products resistant to sterilisation with steam in compliance with the manufacturer's specifications.

- The autoclave must be installed in a suitable and hygienically controlled environment.

- The device must be connected to the electrical mains.

- The electrical system must meet standard IEC 364-1 and the "National rules on the installation of electrical systems on premises suitable for medical use".

- Check the power supply voltage, indicated on the plate, placed on the rear panel of the autoclave, corresponds to that of your mains. If in doubt, contact your retailer.

- The autoclave plug comes with earthing. Ensure you use a plug equipped with an earthing and position the device in such a way it is not difficult to disconnect the plug.

- Using an extension, ensure the type of wire used is suitable for device absorption.

- Before using the device for the first time, clean the room (as described on page 54).

- Do not rest the autoclave on unstable tables or surfaces.

- Do not attempt to carry out maintenance on the autoclave when the power supply is connected; disconnect the power supply cable from the electrical mains, before working on the machine.

- If an error message occurs during a sterilisation cycle, **the cycle must be repeated.**

- To avoid possible burns, the trays must be removed with the specific pincers.

- If the device accidentally falls it may be deformed, which would compromise its good functioning; it should be fully checked by technical support.

- The autoclave works with distilled water (see technical specifications on page 58).

Any other substance used could compromise device functioning. If an inappropriate liquid is used, consult technical support.

- It is forbidden to use any acid to clean any part of the autoclave. In particular, use of hydrochloric acid is forbidden and any detergent containing chlorine. Use of these substances could compromise the mechanical integrity of the device.

- The unit discharges air through a hole in the bottom rear of the device.

- If the appliance is used in a manner not specified by the manufacturer in this manual, the protection provided by the appliance may be impaired.

SAFETY STANDARDS

No harmful substances are used in the sterilization process. We recommend you avoid contact with the content of the discharge tanks.

- Introduce in the sterilization chamber only products resistant to heat that do not release hazardous substances. Always check that the products can be sterilized consulting the technical specifications provided by the manufacturer.
- Never introduce in the chamber substances or products containing heavy metals (lead, mercury, bromine, hexavalent chromium). The release of these products may contaminate the hydraulic circuit of the autoclave and damage it permanently.
- Service and repairs should be carried out by authorized staff of Faro S.p.A. The list of centers authorized to carry out service and repairs is available on the website www.faro.it in the section TECHNICAL SUPPORT/LIST OF TECHNICAL SUPPORT CENTERS. The use must request the intervention of authorized centers only.
- Autoclaves should be used exclusively by staff who undergone suitable training on how to use this manual. Training should be provided by the medical centre in which the autoclave is installed.

- NEVER REUSE THE USED WATER.

- All the objects should be decontaminated and carefully cleaned and dried before being sterilized.
- If there are clustered instruments, such should be divided or however arranged in a position as ventilated as possible.
- Shirts or other reusable fabric should be washed and dried after use and before sterilization to remove any organic material and increase the useful life of the fabric by restoring its natural water content (the level of moisture).

The device, in its packaging, supports the following environmental conditions for a period of no more than 15 weeks:

Environmental temperature from -20°C to +70°C

Relative humidity from 10% to 90%

Atmospheric pressure from 500 to 1060 mBar.

The device must be used in the following environmental conditions:

- indoor use
- altitude up to 3000 m
- temperature from 5° to 40°C
- max relative humidity 85%
- max variation in mains voltage $\pm 10\%$
- installation category (surge category) II
- pollution level I

Do not carry out maintenance operations other than those outlined in the manual. Any intervention not indicated in the manual could compromise the safety aspect, planned for the device. For all operations not specified, contact technical support. N.B. Non-compliant use with that indicated in the user manual could prejudice the safety level of the device. Do not make changes to the device.

The product must be used according to specifications in this manual, in compliance with planned use.

The user is responsible for installation and maintenance and fulfilling the laws relevant to the device. Faro cannot be held liable for damage to property and/or people caused by the incorrect use of the product as described in this manual. Non-compliance with the specifications in this manual means any manufacturer liability and product guarantee becomes null and void.

CHARACTERISTICS

THE DEVICE IS INTENDED FOR STERILISATION USING WET SATURATED STEAM.

Use the vacuum pump to obtain an optimal sterilisation process with correct interaction between the steam and the load processed also in the cavities.

STERILISATION PROGRAMS

TYPE B

- 121°C or 134°C – **cycle type B standard**
- 134°C – **cycle type B Prion**
- 134°C B Flash – **cycle type B fast**

TYPE S (optional)

- 121°C – **porous materials**
- 134°C – **wrapped material**
- 134°C – **unwrapped material**

TEST PROGRAMS

TYPE B

- **Vacuum test:** checks the capacity to maintain the vacuum
- **Bowie & Dick test:** checks steam penetration in the porous load
- **121 Helix test**
- **134 Helix test:** check steam penetration in hollow loads

TYPE S

- **Vacuum test**
- **Bowie & Dick test**

The autoclave is available in 2 capacity variants in the sterilisation chamber:

- 16
- 21

The autoclave is available with the following data saving or data print-out devices:

- external printer (data print-out only)
- built-in printer (data print-out only)
- memory card (saves data with print-out via dedicated sw).

AUTOMATIC LAST ADVANCED CYCLE FUNCTION

When the cycle is complete, with the door closed, the autoclave goes in stand-by after 30 minutes. When the temperature of the chamber reaches 60° C the drying phase is automatically repeated.

WATER FILLING SYSTEMS

Automatic: with the help of a pump inside the autoclave.

Mains connection: possibility of connecting a demineralisation system, to increase autonomy. Installation must be carried out by an authorised technician.

Manual: alternatively, the device is equipped with emergency manual filling on the rear.

CHARACTERISTICS

AUTOMATIC CONTROL OF STERILISATION PARAMETERS

Pressure, temperature and time parameters are continuously monitored during the sterilisation phase using an automatic process control system.

BIO-X FILTER

Releases sterile air into the chamber during the drying phase and pressure balancing phase.

USB PORT

Unpowered port to be used exclusively by qualified engineers. Use a standard USB cable with connectors type A - type B with max length 3m.

SERIAL PORT - RS232

Port dedicated exclusively to FARO printer connection. Use the dedicated serial cable supplied with the FARO printer.

OLED display

A menu system enables complete management of the autoclave.

DESCRIPTION OF PARTS

SEE FIG. 1

1. Switch
2. Knob to open – close door
3. LCD display
4. Button (\\) - Reset
5. Button (←) - Enter
6. Button (-) – Scroll up
7. Button (+) – Scroll down
8. USB connector
9. Bio-X filter
10. Built-in printer compartment/memory card slot/external printer compartment

SEE FIG. 2 – FRONT VIEW

11. Coupling for discharge tank discharge
12. Coupling for filling tank filling
18. Handles

SEE FIG. 3 – REAR VIEW

13. Coupling for filling tank discharge in case of service
14. Coupling for manual filling of filling tank
15. Water filling joint from the demineralizer to connect to the water mains
Pressured permitted from 0.2 : 2.5 bar
16. Coupling to discharge the discharge tank

SEE FIG. 4

17. Door adjustment button
19. Porthole

TYPE TEST ON AUTOCLAVE MODEL SK07

Type test conducted according to standard EN13060 Type B cycles

Pressure dynamic in sterilisation chamber
Air leakage
Vacuum chamber
Solid load
Small porous load
Complete porous load
Hollow load A
Multiple wrapping
Drying, solid load
Drying, porous load

Type test carried out according to standard EN13060 Type S cycles

	121°C Porous Load	134°C Wrapped	134°C Unwrapped
Air leakage	X	X	X
Pressure dynamic in sterilisation chamber	X	X	X
Vacuum chamber	X	X	X
Solid Load unwrapped			X
Solid Load individual wrapping		X	
Small Porous Load individual wrap	X		
Drying Solid Load unwrapped			X
Drying Solid Load individual wrap		X	
Drying Small Porous Load individual wrapping	X		

INSTALLATION

THE DEVICE MUST BE INSTALLED BY SPECIALIST TECHNICIANS.

Check the installation space ensures the minimum space indicated can be maintained.


FIRST SWITCH-ON

Qualify the installation according to the following test table:


- Vacuum test
- Automatic control test
- Verification of the measuring instruments of the sterilizer
- Overheat cutout test
- Thermometric test for small load (*)
- Small load dryness test (*)
- Steam penetration test for porous loads (*)
- Thermometric test for solid load (*)
- Solid load dryness test (*)
- Steam penetration test for solid loads (*)
- Thermometric tests (load selected on user requirements)
- Microbiological tests (load selected on user requirements)

(*) To be carried out with cycles of 121°C and 134°C

MODEL SIZES 16-21

The door has the following open dimensions:

The bottom of the autoclave packaging has two gripping points.

The short movement operation of the device must be carried out by two people; the gripping points to lift the autoclave (handles) are located on the top of the device (18 Fig. 2). For prolonged movement of the device you must insert the pallet supplied under the device.

- Take the autoclave out of its packaging.
- Position the autoclave on a surface (table or unit) using the specific handles to help you.
- Remove the handles and insert the cover plugs. Keep the handles and the relevant screws for future movement.
- Open the door, turning knob 2 in an anti-clockwise direction (Fig. 1) and remove the accessories supplied from the chamber.

Check the supply includes:

- Documentation (user manual)
- The tray support
- 3 trays
- 1 tray extraction pincers
- 1 filling tube
- 1 discharge tube
- 1 funnel
- 4 plugs

For versions equipped with Memory Card:

- Memory card inserted in dedicated slot
- USB memory card reader for PC
- CD for reading and printing data from memory card.

Use a level to check the support base is level in the two directions, otherwise level the autoclave using the front adjustment feet.

Check the power supply socket is earthed.

Check the power supply voltage corresponds to that on the autoclave plate.

Insert the autoclave plug while paying attention that the wire is not crushed and that the power supply plug is easily removed.

If the water is not filled from external canisters, the autoclave MUST NOT BE CONNECTED directly to the water mains, but to a demineralizer that permits the required water characteristics to be obtained (see page 62). The demineralizer must prevent the backflow of water to the mains.

INSTRUCTIONS FOR USE

To use the autoclave, proceed as follows:

- Open the door by turning the knob (2 Fig.1) in an anti-clockwise direction.
- Press the main light switch (1 Fig. 1).

USING THE MENUS


The main menu is composed of 3 sub-menus:


- 1 - Programs
- 2 - Tests
- 3 - Options





To access the various sub-menus, move the selection indicator by pressing the (+) or (-) keys and then confirm your choice with the  key.

Press the (\\) key to cancel the operation.

MAIN MENU

FILLING THE FILLING TANK

- Only use distilled water (see technical specifications on page 62).
- When the LCD display shows the message “Proceed to fill H2O” fill the filling tank with distilled water, choosing from one of two of the following modes:

AUTOMATIC FILLING

- Insert the tube part, with the quick coupling, in the specific water filling joint (**12 Fig. 2**) and insert the other end in the canister.
- Press the (+) button to start the water filling pump.
- Press the (-) key (\\) to stop the operation.
- Filling is automatically blocked once the maximum level of the filling tank is reached.
- **Automatic filling of feeding tank from water mains.**

MANUAL FILLING

- Insert the tube part with the quick coupling in the hole (**14 Fig. 3**).
- Insert the funnel into the other end of the tube and pour max 2 litres of water.

EMPTYING THE DISCHARGE TANK

- Insert the side of the tube without the quick coupling in a container to collect the discharged water.
- Insert the quick coupling in the hole **11 Fig. 2** and position the container under the autoclave to empty the tank.
- **Water must be discharged in compliance with the standards in force on waste disposal.**
- **Automatic emptying of discharge tank (16 fig. 3)**

PLACING MATERIAL TO STERILISE IN THE CHAMBER

Wrap the load very carefully determining the compatibility of the medical devices to wrap with the pre-chosen wrapping.

The material used for packaging must ensure:

- **Compatibility with the sterilisation process.**
- **Compatibility with the labelling system.**

Always insert a chemical indicator with the load to verify the positive outcome of the sterilisation thermodynamic cycle.

Arrange the well-cleaned material on the trays, distributing it uniformly.

The chamber can have an inserted maximum of:

	Model 16 L	Model 21 L
wrapped material	2,5 kg	3,5 kg
unwrapped material	4 kg	6 kg
porous load	1,3 kg	1,6 kg

The wrapped objects must be placed on trays with the transparent part of the wrapping turned downwards.

“PROGRAMMED START” FUNCTION

From the menu the user can access the sub-menus **“Instant start - Programmed Start”**

- a-** Select **Instant start** to access the sterilization cycle or the test cycle.
- b-** Select **Programmed Start** to access the menu for setting the start time. The user can set the hours and the minutes by entering every digit using the keys +/- selecting with the Enter key . To confirm the setting go to **Save** and confirm with the Reset key.

Once you have set the Start time, access the menu for selecting the sterilisation cycle or the test cycle.

STERILISATION PROGRAM SELECTION

The pre-heating time, with the machine cold, is about 10 minutes. By selecting the “**Programs**” sub-menu you access the list of sterilisation cycles that the autoclave can execute:

CYCLE B

- Cycle 134°C Standard
- Cycle 134°C Prion
- Cycle 121°C Standard
- Cycle 134°C B Flash*

CYCLE S

- Cycle 121°C porous material
- Cycle 134°C wrapped
- Cycle 134°C unwrapped

To scroll through the cycles and move the selection indicator press the (+) or (-) keys and then confirm your choice with the **OK** key. Press the “\” key to cancel the operation.

STANDARD CYCLES

CYCLE B

- 134°C: 4 min. sterilisation plus 14 min. drying
- 134°C B Flash*: 3,5 min. sterilisation plus 8 min. drying
- 121°C: 16 min. sterilisation plus 14 min. drying

CYCLE S

- 121°C: 16 min. sterilisation plus 8 min. drying
- 134°C wrapped: 4 min. sterilisation plus 8 min. drying
- 134°C B unwrapped: 4 min. sterilisation plus 6 min. drying

* maximum recommended load per version 16: 0.5 Kg wrapped load – 0.2 kg porous load. For version 21: 0.7 Kg wrapped load – 0.3 kg porous load.

NOTE: For each sterilisation cycle available, the manufacturer must carry out Works Tests in compliance with EN 13060 (par. 7.3.1). The cycles available for which the manufacturer has not carried out Works Tests must undergo Installation Tests according to EN 13060 par. 7.4.

PRION CYCLE

- **Ciclo at 134°C: 18 minutes sterilisation plus 14 minutes drying.**

HOW TO USE THE MEMORY CARD

To use the Compact Flash card, comply with the instructions supplied on the CD. To download the data, use the specific program supplied.

“LAST ADVANCED CYCLE” FUNCTION

The “**LAST ADVANCED CYCLE**” function is activated automatically for all the programs. 30’ after the end of the sterilisation cycle, the heating elements switch off and a **second drying cycle automatically takes place** when the temperature inside the chamber reaches 60°C.

ISTERILISATION TIMES INTERVENTION

You can increase the sterilisation times; the operation must be carried out by a specialist technician.

“ALERT DISPLAY” FUNCTION

Quick flashing: from max to min every second when there are 4 minutes to cycle completion.

“Stand -by display” FUNCTION

Slow flashing: from max to min every 6 seconds during stand-by (up to 30 minutes from cycle completion).

“Save page on display” FUNCTION

FARO logo slider: after stand-by.

“PRESSURE UNIT OF MEASUREMENT” FUNCTION

Use the Options-Setup-Special-Pressure UM menu to select the Pressure unit of measurement that you want to be displayed.

You can choose from:

- KPa (Kilo Pascal)
- Bar (Absolut)
- Bar (relative)

“DISPLAY CONTRAST” FUNCTION

Use the Options-Setup-Special-Display Contrast menu to select the desired intensity of the display.

You can choose from:




- maximum
- minimum
- default

During setup is displayed a number that varies from 50 to 255 that indicates the level of contrast.





POSSIBILITY TO ENABLE/DISABLE THE FLASHING ALERT AND STAND-BY FUNCTIONS

Use the Options-Setup-Special-Display Contrast menu to enable/disable the flashing function on display.




CYCLE "B" DIAGRAMS



CYCLE "B" DIAGRAMS

CYCLE "S" DIAGRAMS

CYCLE "S" DIAGRAMS

EXAMPLE OF PRINTED VALUES

PRINTING SEQUENCE:

Logo Faro or FARO S.p.A.

Ornago (Mi)

SNXX STXXXX ← Autoclave serial no.

Program: ← type of cycle or type of test

Cycle nr XXXX ← no. of cycle

Start: ← date and time cycle starts

DD/MM/YY HH:MM:SS

PRESS TEMP TIME ← unit of measurement
(kPa) (°C) (H:MM:SS) (time calculated from start of cycle)

Fractional Vacuum

XXXX XXXX X:XX:XX ← Pressure, temperature and time of 1st vacuum phase
+++ ← Top-ups during pressure increase

XXXX XXXX X:XX:XX ← Pressure, temperature and time of 2nd vacuum phase
+++ ← Top-ups during pressure increase

XXXX XXXX X:XX:XX ← Pressure, temperature and time of 3rd vacuum phase
+++ ← Top-ups during pressure increase

Sterilization

XXXX XXXX X:XX:XX ← Pressure, temperature and time on starting sterilization

XXXX XXXX X:XX:XX ← Pressure, temperature and time during sterilization

XXXX XXXX X:XX:XX ← Pressure, temperature and time during sterilization

XXXX XXXX X:XX:XX ← Pressure, temperature and time during sterilization

XXXX XXXX X:XX:XX ← Pressure, temperature and time at end of sterilization

+++ ← Top-ups during sterilization

P_{MAX}=XXXXkPa T_{MAX}=XXXX°C ← Maximum pressure and temperature during sterilization

P_{MIN}=XXXXkPa T_{MIN}=XXXX°C ← Minimum pressure and temperature during sterilization

Discharge

X:XX:XX ← Discharge phase start time

Dry

X:XX:XX ← Drying phase start time

Press. Balance

X:XX:XX ← Pressure balancing phase start time

End:

X:XX:XX ← Cycle total time

dd/mm/yy hh:mm:ss ← Date and time cycle ends

OK

Operator signature

TEST PROGRAMS SELECTION

By selecting the “test” sub-menu you access the list of test cycles that autoclave must execute:

- 1 - Vacuum Test
- 2 - 134 Helix Test
- 3 - 121 Helix Test
- 4 - Bowie&Dick Test
- 5 - Exit

The selection indicator is positioned on “Vacuum Test”. Select the desired program by moving with the (+) and (-) keys and confirm with .

VACUUM TEST (Leakage test)

To execute this function the chamber must be cold and contain no instruments.

BOWIE & DICK TEST

There must be no tools in the chamber to carry out this test. The Bowie & Dick “package” to test the porous load must be positioned on the lower tray.

HELIX TEST

There must be no tools in the chamber to carry out this test. The control tool to test the “Hollow A” load must be positioned on the lower tray.

OPTIONS SELECTION:

By selecting the “Options” sub-menu you access the list of configuration sub-menus for the machine:

- 1 - Configuration
- 2 - Fill water
- 3 - Print cycle
- 4 - Delete memory
- 5 - Exit

Configuration sub-menu:

By selecting the “Configuration” sub-menu you access the list of configuration options for the machine:

- 1 - Basic
- 2 - Advanced
- 3 - Special
- 4 - Exit

Basic selection:

By selecting the “Basic” sub-menu you access the list of configurations for the machine:

- 1 - Language
- 2 - Date
- 3 - Time
- 4 - Exit

Advanced selection:

Only an authorised technician can access this menu. By selecting the “Advanced” sub-menu you access the list of configurations for the machine:

- 1 - Fill water options
- 2 - Saving
- 3 - Exit

Special selection:

Only an authorised technician can access this menu. By selecting the “Special” sub-menu you access the list of configurations for the machine:

- 1 - Environmental Pressure
- 2 - LCD Contrast
- 3 - Password Management
- 4 - Exit

MAINTENANCE

Operations	Frequency			
	Daily 5 cycles	Monthly 50 cycles	Quarterly 500 cycles	Annual 1000 cycles
Porthole cleaning	X			
Porthole gasket cleaning	X			
Chamber and accessory cleaning		X		
Front panel cleaning		X		
Chamber filter cleaning		X		
Porthole door adjustment				
Bacterial filter replacement			X	
Porthole gasket replacement				X

PERIODS CONTROLS

Operations	Frequency			
	Before first use	Daily	Annually	Every two years
Vacuum Test EN 17665		X		
Bowie & Dick Test 17665		X		
Helix Test EN 17665		X		
Front panel cleaning - Safety earthing continuity - Insulation test - Casing dispersion current				X
Initial validation EN 17665	X			
Periodic validation EN 17665			X	
Safety valve check			X	

! Before proceeding with the operations indicated below, always switch off the main switch (1 Fig. 1) and check the chamber cools.

Use is forbidden of hydrochloric acid and detergents containing chlorine to clean any part of the autoclave.

- Cleaning the sterilisation chamber and the accessories inside

Check: Monthly or every 50 cycles.

Use a soft cloth soaked in water so that the protective surface treatment of the chamber is not damaged. Rinse with plenty of water and, at the end of the operation, dry. If necessary, dampen the cloth with white vinegar to remove any limescale.

You are however advised to clean the chamber only if strictly necessary.

If limescale deposits are present, immediately check what is causing them. The formation of white patches at the bottom of the chamber means poor quality water is being used or water with incorrect characteristics compared to those required.

- Porthole cleaning

Check: Daily or every 5 cycles.

Use a soft cloth soaked in water to avoid damaging the protective treatment of the porthole.

Rinse with plenty of water and when finished, dry. During cleaning DO NOT USE metal type or particularly abrasive cloths. Use of the latter could cause damage to the outer surface of the porthole.

- External front panel cleaning

Check: Monthly or every 50 cycles.

Use a soft cloth soaked in alcohol, taking care it does not penetrate the inner parts of the autoclave. It is strictly forbidden to use detergents with a Trilene, benzene or turpentine base or solvents of any kind.

- Porthole gasket cleaning

Check: Daily or every 5 cycles.

Clean the gasket with a soft cloth dampened with water or white vinegar, to remove any limescale. Rinse well if you have used vinegar.

This procedure must be carried out to remove impurities that can cause pressure loss in the sterilisation chamber and possible gasket breakage.

- Chamber filter cleaning

Check: Monthly or every 50 cycles.

- Remove the tray supports with trays from the chamber
- Turn the filter in the chamber and pull towards yourself to remove it.
- Place the chamber filter in water or vinegar, wait a few minutes and then with a small brush carefully clean the filter area.
- Rinse the filter unit with plenty of water.
- To remove the filter, proceed with the operations in the opposite order.
- Re-position the filter on the back of the chamber.

If the cleaning operation is not satisfactory, call Technical Support.

- Porthole door adjustment

Check: When necessary. Normally no earlier than six months from the porthole gasket change.

The autoclave closure system requires occasional adjustment, due to normal adjustment of the mechanical parts and wear on the gasket.

Imperfect hold on the gasket could prevent correct execution of the sterilisation cycles and therefore compromise the good outcome.

Always work with the machine off and cold.

Operation:

- Keep button 17 Fig. 4 pressed
- Turn porthole 19 Fig. 4 in an anti-clockwise direction for a quarter of a turn
- Release the button and continue rotating porthole 19 Fig. 4 until it locks
- Check the door closes regularly. If the handle is too difficult to close, turn slightly from the opposite side while keeping button 17 Fig. 4 pressed.
- Run a test cycle to check adjustment was efficient.

This procedure usually should take place on the first slight leak or if the gasket has difficulty holding. Before adjusting the porthole closure system, clean the porthole gasket as indicated above and repeat the sterilisation cycle.

- Bacterial or Bio X filter replacement

Check: Quarterly or every 500 cycles (replacement is advised)

Operation:

- Open the door on the printer/CF Card compartment on the autoclave.
- Unscrew the filter (9 Fig. 1) by turning in an anti-clockwise direction.
- Screw the new filter in, in a clockwise direction, paying attention to completely screw it in.

- Porthole gasket replacement

 **Always work with the machine off and cold.**

Check: Annually or every 1000 cycles.

Operation:

- Open the door of the steriliser
- Manually remove the lip gasket, grabbing the lip with two fingers (pay attention as to how the gasket lip is assembled).
- Clean the gasket compartment using a cloth dampened with alcohol or a cotton bud.
- Dampen the new gasket with water or alcohol.
- Insert the new gasket in the compartment with the lip turned outwards; distribute the gasket evenly on the circumference and with your finger evenly press on the entire perimeter of the gasket.
- When the procedure is complete, check there are no badly inserted points or the gasket is not protruding. If the gasket is protruding, remove the gasket and insert again.

 **Correctly position the gasket to avoid leaks. The gasket must be inserted over the entire circumference without bulges or deformations.**

ANNUAL MAINTENANCE

After 1000 cycles or a year after installing the autoclave.

On reaching cycle no. 1000 a memo appears to carry out maintenance. Maintenance includes: replacement, cleaning and checking the parts, as indicated in the relevant table. Annual maintenance can only be carried out by a FARO qualified technician and a member of the FARO Technical Support Service network.

Intervention includes the following:

Replacing:

- Porthole gasket
- Bacterial or Bio X filter
- Filter in chamber
- Water filter
- Replace vacuum pump membrane kit (if necessary)

Check:

- Check the door safety closure system
- Check the porthole door adjustment
- Check and clean the water filling and discharge tank
- Check the hydraulic connections
- Check the electrical connections
- Check calibrated pressure and temperature with SW
- Check the vacuum pump is working properly.

ANNUAL MAINTENANCE


Check:

- CSafety valve
- Electrical dispersion

When the qualified technician has finished, a certificate is issued to state annual maintenance of the autoclave took place.

To ensure the guarantee of the device is valid, comply with the intervals specified or indicated in this manual and in the additional sheets and only use ORIGINAL SPARE PARTS.

DISPLAY SYMBOLS

CLOCK SETTINGS

Date and Time Adjustment:

From the main menu select “Options”.
 From the options menu select “Configuration”.
 From the configuration menu select “Basic”.

To change the date:

From the Basic menu select “Date”:

Press the “+” key
to increase
selection.

Press the key “-”
to decrease
selection.

Press the “Enter” key”
to move with the selection.

Select “Save” and press
the “Reset” key
to save the settings.

To change the time:

From the Basic menu select “Time”:

Press the “+” key
to increase
selection.

Press the “-” key
to decrease
selection.

Press the “Enter” key
to move with the selection.

Select “Save” and press
the “Reset” key
to save the settings.

CHANGE PRINTING PAPER ROLL

To change the paper roll, proceed as follows:

1. Open the printer lid by pressing the OPEN key as indicated (see **fig. A**).
2. Position the paper roll inside the chamber, respecting the indicated rotation direction (see **fig. B**).
3. Pull the paper until it exits the chamber and close the door (see **fig. C**). The door locks.
4. Tear away the excess paper using the jagged edge.

fig. A

fig. B



fig. C

DOOR RELEASE PROCEDURES

MACHINE RELEASE PROCEDURE IF ERROR OR ALARM MESSAGE APPEARS

If an Error or Alarm message appears, the lower part of the display will show the **“Release machine”** message.

To restore the autoclave **SIMULTANEOUSLY PRESS** the keys:


MACHINE RELEASE PROCEDURE IF “OPEN DOOR” MESSAGE APPEARS WITH SUCCESSFUL CYCLE

If the door is locked, when the “Open door” message appears, do as follows:

- 1) Bring the opening handle to 90° as in the figure and keep it in this position until the procedure is complete.


- 2) Press and keep the key below pressed for at least 2 sec.

- 3) The autoclave display shows:
AL.015 - Reset Pressed

- 4) Wait until the lower part of the display shows the message **“Release machine”**.

To restore the autoclave **SIMULTANEOUSLY PRESS** the following keys:

ERRORS DURING CYCLES

COD.	MESSAGE	DESCRIPTION	ACTION
ER001	Safety Cut Off	safe temperature exceeded in chamber or band	Contact FARO technical support, stating the data on the autoclave display.
ER002	Top-ups	maximum number exceeded of permitted top-ups	Check the load quantity put in the autoclave.
ER003	Pressure balancing	pressure balancing not reached	Replace the Bio-X filter.
ER004	No Ready	ready temperature not reached	Contact Faro technical support.
ER005	Timeout 1	final pressure not reached	Check the load quantity put in the autoclave.
ER006	Vacuum	second vacuum threshold not reached	Check the porthole status and door gasket status.
ER007	Timeout 2	sterilisation conditions not reached	Check the load quantity put in the autoclave and the porthole and door gasket status.
ER008	T chamber High	chamber temperature over permitted limit	Check chamber probe position.
ER009	T chamber Low	chamber temperature under permitted limit	Check chamber probe position.
ER010	High Pressure	pressure over permitted limit	Contact Faro technical support.
ER011	Low Pressure	pressure under permitted limit	Check the load quantity put in the autoclave.
ER012	T theoretic High	theoric temperature over permitted limit	Check the position of the chamber probe. If necessary, contact Faro technical support.
ER013	T theoretic Low	theoric temperature under permitted limit	Check the chamber probe position.
ER014	Delta Temperature	difference between chamber and theoric temperature $> +/ - 2K$	Check the chamber probe position.
ER015	Vacuum Test: vacuum	vacuum temperature not reached during vacuum test	Check the porthole and door gasket status.
ER016	EN 13060	non-compliance with parameters in 300 sec	Check the porthole and door gasket status.
ER017	EN 13060	non-compliance with parameters in last 600 sec	Check the porthole and door gasket status.
ER018	EN 13060	non-compliance with temperature limits	Run Vacuum Test on cold machine.
ER019	Door open	no door closure	Check the door has closed.
ER020	Water Fill, Timeout	filling tank not filled	Fill the filling tank.
ER021	Empty Discharge Tank	filling tank not emptied	Empty the discharge tank.

To reset the fault and restore the system, do as follows:
 Simultaneously press the “ ← ”, “ + ” and “ - ” keys and wait for the procedure to end.

GUIDE TO PROBLEMS

MACHINE FAULTS

COD.	MESSAGE	DESCRIPTION	ACTION
AL001	Chamber Probe Fault	chamber probe open	Contact Faro technical support.
AL002	Band Probe Fault	chamber probe in short circuit	Contact Faro technical support.
AL003	Guasto Sonda Fascia	band probe open	Contact Faro technical support.
AL004	Guasto Sonda Fascia	band probe in short circuit	Contact Faro technical support.
AL005	Water Fill. Timeout	maximum time exceeded to fill water	Check filling tube and presence of water in can.
AL006	Pressure Transducer	pressure transducer fault	Contact Faro technical support.
AL007	Door Micro	unplanned door micro contact opening	Contact Faro technical support.
AL008	Safety Micro	unplanned door safety micro contact opening	Contact Faro technical support.
AL009	Blackout	blackout	Check if shut-off due to electrical system.
AL010	Safety Pressure	pressure over safe limit	Contact Faro technical support.
AL011	Corrupt Memory	corrupt memory	Contact Faro technical support.
AL012	Safety Micro	No door safety inserted	Contact Faro technical support.
AL013	PUSH Error	Software error	Contact Faro technical support.
AL014	POP Error	Software error	Contact Faro technical support.
AL015	Reset Pressed	Reset key pressed	Proceed with autoclave release procedure.
AL016	Safety Micro	No door safety disconnection	Contact Faro technical support.

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Protection rating against penetration of solid and liquid particles	: IP20
Power cord features	: Operating temperature 60°C
Fuses replacement	: The operation should be carried out only by qualified engineers

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GUARANTEE CERTIFICATE

GB

Faro offers the final customer a **24 month** guarantee starting from the date of purchase.

Repairs under guarantee must be performed at FARO; expenses and transport risks are at the risk of the purchaser.

Repair under guarantee is considered valid only when:

- **all sections of the certificate have been filled in and sent in advance to FARO by Fax (039.6010540).**

The guarantee covers faults due to the bad quality of the material or manufacturing defects; in the case of valid claims, the guarantee covers free repair or replacement. **Claims for damages and/or interest are excluded.** The guarantee is not considered valid, at the sole discretion of FARO, if the fault is due to tampering, damage, incorrect use, improper maintenance and normal wear and tear.

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FARO SpA si riserva il diritto di modificare, senza preavviso, le caratteristiche indicate nel presente manuale.
FARO SpA reserves the right to change the specifications of this equipment without notice.
FARO SpA se reserve le droit de modifier, sans pr avis, les caract ristiques dans ce manuel.
FARO SpA beh lt sich rechtvor, jederzeit stillschweigend technische oder bauliche nderung vorzunehmen.
FARO SpA se reserva el derecho de modificar sin aviso previo la caracteristicas incluidas en el presente manual de uso.